


Non-Adhesive Ostomy Systems



Product Catalog

Table of Contents

Urostomy Systems _____	2 - 4
Ileostomy Systems _____	4 - 9
Colostomy Systems _____	9 - 10
O Rings _____	11
Belts _____	12
Pouch Comfort Panel _____	13

Nu-Hope Laboratories, Inc.

The history of **Nu-Hope** is an interesting one of tragedy turned triumph because Nu-Hope was born out of necessity. The founder, Edmund Galindo, suffered from chronic kidney disease. Unaware of the severity of his condition, he volunteered for the war effort in 1941. Doctors there were so alarmed they told him he only had one more year to live! He survived that prediction by 18 years. Eventually medical treatment was no longer effective and a urinary diversion was performed in January of 1958. Back then there were no ET's to turn to for help and there were many frustrating months of leakage, sore skin, and constant laundry. However, Edmund was not about to give up. Being of keen mind, he experimented with new appliance designs and by the middle of 1958 he succeeded. Ed proved the worthiness of his new appliance system on himself, returning to his active life of business, golf and swimming. The doctors in the Los Angeles area were so impressed with his success that they urged him to manufacture and market his appliances so that others might benefit.

At first, supplies and appliances were freely handed out. Soon it became apparent that in order to continue helping others a company had to be formed. Thus, in 1959, Nu-Hope Laboratories, was founded. The name was derived from Mrs. Galindo's first name "Hope", and the feeling that they were offering "new hope" to others in need. Unhappily, Edmund Galindo's health continued to deteriorate and on November 11, 1959, he passed away. Determined to carry on the founder's legacy, Hope Galindo and 19 year old son Eugene continued with the business. In 1974 Edmunds daughter, Louise, became an active member of the family business for 20 years.



During the pioneering years, Gene Galindo visited and fitted hundreds of ostomy patients in hospitals and homes and learned much about stomas and stoma care. Dealing with the many different shapes and sizes presented new challenges. Stoma placement, skin conditions, and the new ostomate's acceptance of their ostomy were all issues that needed to be addressed. Gene continued in his father's footsteps by developing many news ostomy appliances and related accesories. Enterostomal Therapists did not exist at that time and Nu-Hope partially filled that need, working closely with the medical profession and ostomy associates.

Currently there are over 4,000 E.T.'s and several Wound Ostomy Continence Nurse accredited Educational Programs. Nu-Hope continues to work closely with E.T.s/WOC nurses to provide innovative products and services for their special needs. Over the years Nu-Hope has built a reputation of manufacturing innovative, high quailty, reliable ostomy appliances and accessories.

Ordering Information

All orders should be directed to the Customer Service Department, by mail, telephone, fax, or email. Toll free line 800.899.5017 is available weekdays (excluding national holidays) between the hours of 8:00 a.m. and 4:30 p.m. PST. Fax number: 818.899.2079. Send email inquiries to: info@nu-hope.com.

Shipment and Terms

All stock products are shipped UPS Ground Service unless otherwise specified at the time of ordering. Terms are F.O.B. origin, prepaid and charged. Customers may request priority shipping. Orders shipping outside the USA will be charged full shipping charges, shipping charges to be either prepaid and added to the invoice or shipped freight collect. Any special requests may be subject to a minimal \$5.00 service charge added to the invoice.

Product Warranty

We warrant that, at the time of manufacture, this product was prepared and tested and meet our specifications.

Because of biological differences in individuals, no product is 100% effective under all circumstances. In addition, because we have no control of the condition under which the product is used, diagnosis of the patient, the method of use or administration, and handling of the product after it leaves our possession, we do not warrant either a good effect or against an ill effect following the products use.

The foregoing warranty is exclusive and in lieu of all other warranties either written, oral or implied (including any warranties of merchantability or fitness for purpose). No representative of the company may change any of the foregoing and the buyer accepts the product subject to all terms herein.

Returned Goods Policy

Authorization is required before returning merchandise. Items returned without a Return Authorization number will be returned freight collect at our option. For all returns of any stock item, contact Nu-Hope by calling 800.899.5017. Requests must specify reason for return (e.g. ordering error, shipping error), catalog product number and quantity and invoice number with date of invoice and/or purchase order number. Returned products must be accompanied by original order invoice, packing slip or facsimile. Returned items are subject to a reprocessing charge up to 20% of original item cost. Full credit will be given to the customer for an item returned due to a processing error by Nu-Hope. No credit will be given for an item returned unsuitable for restocking. Returned products must be accompanied by an RA number, be undamaged, properly packaged for shipping, and should not have been used.

Due to variations within product lines, Nu-Hope reserves the right to designate stock items.

Product Changes

Nu-Hope reserves the right to change or discontinue any product without notice. If the product has been discontinued, we can recommend alternatives.

Product Storage


Nu-Hope Non-Adhesive products should be stored in a dark, cool and dry place. Extended periods of exposure to light have resulted in degradation of some vinyl materials and should be avoided. Extended periods of elevated temperature can reduce a product's shelf-life and should also be avoided.

New Customers

Customers who order an ostomy product for the first time and discover that it is not the correct size or does not meet their need, will have 30 days from the date the product is delivered, to receive a full refund of the cost of the product.

Urostomy Systems

Used to provide urostomates a reusable system for urine collection and drainage. Supplied non-sterile. Intended for single patient use only.



Urostomy Systems

Right Stoma Location

Product Number	Description	O-Ring Size	Pouch Size
EV5000-004	Sets consist of 1 (one) pouch, one o-ring and one belt	Small ($\frac{7}{8}$ ")	X-Small (300ml)
EV5000-000			Small (500ml)
EV5000-001			Medium (800ml)
EV5000-002			Large (1100ml)
EV5000-4S1	Sets consist of 1 (one) pouch, one o-ring and one belt	Small ($\frac{7}{8}$ ") Extra Tall	X-Small (300ml)
EV5000-0S1			Small (500ml)
EV5000-1S1			Medium (800ml)
EV5000-2S1			Large (1100ml)
EV5010-004	Sets consist of 1 (one) pouch, one o-ring and one belt	Medium ($\frac{7}{8}$ " - $1\frac{1}{8}$ ")	X-Small (300ml)
EV5010-000			Small (500ml)
EV5010-001			Medium (800ml)
EV5010-002			Large (1100ml)
EV5010-0S1	Sets consist of 1 (one) pouch, one o-ring and one belt	Medium ($\frac{7}{8}$ " - $1\frac{1}{8}$ ") Extra Tall	Small (500ml)
EV5010-1S1			Medium (800ml)
EV5010-2S1			Large (1100ml)
EV5010-0S2	Sets consist of 1 (one) pouch, one o-ring and one belt - Special Belt Used	Medium ($\frac{7}{8}$ " - $1\frac{1}{8}$ ")	Small (500ml)
EV5010-0S6			Medium (800ml)
EV5010-1S6	Metal Adjusters Used	Medium ($\frac{7}{8}$ " - $1\frac{1}{8}$ ")	Medium (800ml)
EV5020-004	Sets consist of 1 (one) pouch, one o-ring and one belt	Large ($1\frac{1}{8}$ " - $1\frac{3}{8}$ ")	X-Small (300ml)
EV5020-000			Small (500ml)
EV5020-001			Medium (800ml)
EV5020-002			Large (1100ml)
EV5020-4S6	Sets consist of 1 (one) pouch, one o-ring and one belt	Large ($1\frac{1}{8}$ " - $1\frac{3}{8}$ ") Extra Tall	X-Small (300ml)
EV5020-0S1			Small (500ml)
EV5020-1S1			Medium (800ml)
EV5020-2S1			Large (1100ml)
EV5020-0S4	Sets consist of 1 (one) pouch, one o-ring and one belt	Extra Large ($1\frac{1}{4}$ ")	Small (500ml)
EV5020-1S4			Medium (800ml)
EV5020-2S4			Large (1100ml)
EV5020-2S6	Sets consist of 1 (one) pouch, one o-ring and one belt - Metal Adjusters Used	Large ($1\frac{1}{8}$ " - $1\frac{3}{8}$ ")	Large (1100ml)


Urostomy Systems

Left Stoma Location


Product Number	Description	O-Ring Size	Pouch Size
EV5002-000	SETS consist of 1 (one) pouch, one o-ring and one belt	Small ($\frac{7}{8}$ "	Small (500ml)
EV5002-001			Medium (800ml)
EV5002-002			Large (1100ml)
EV5012-004	SETS consist of 1 (one) pouch, one o-ring and one belt	Medium ($\frac{7}{8}$ " - $1\frac{1}{8}$ "	X-Small (300ml)
EV5012-000			Small (500ml)
EV5012-001			Medium (800ml)
EV5012-002			Large (1100ml)
EV5022-004	SETS consist of 1 (one) pouch, one o-ring and one belt	Large ($1\frac{1}{8}$ " - $1\frac{3}{8}$ "	X-Small (300ml)
EV5022-000			Small (500ml)
EV5022-001			Medium (800ml)
EV5022-002			Large (1100ml)
EV5022-4S1	SETS consist of 1 (one) pouch, one o-ring and one belt	Large ($1\frac{1}{8}$ " - $1\frac{3}{8}$ " Extra Tall	X-Small (300ml)
EV5022-1S1			Medium (800ml)
EV5022-1Z1	SETS consist of 1 (one) pouch, one o-ring and one belt - Special Belt Used	Large ($1\frac{1}{8}$ " - $1\frac{3}{8}$ "	Medium (800ml)

Ileostomy Systems

Used to provide ileostomates a reusable system for fecal collection and drainage.
 Supplied non-sterile.
 Intended for single patient use only.



Emptying Pouch



Ileostomy Systems

STARTER Sets - Right Stoma Location

Product Number	Description	O-Ring Size	Pouch Size
EV6100-004	SET consist of 1 (one) pouch, one o-ring, one belt and one drain tail closure.	Small ($\frac{7}{8}$ "	X-Small (300ml)
EV6100-000			Small (500ml)
EV6110-004		Medium ($\frac{7}{8}$ " - $1\frac{1}{8}$ "	X-Small (300ml)
EV6110-000			Small (500ml)
EV6120-004		Large ($1\frac{1}{8}$ " - $1\frac{3}{8}$ "	X-Small (300ml)
EV6120-000			Small (500ml)

STARTER Sets - Right Stoma Location - *Extra Tall*

Product Number	Description	O-Ring Size	Pouch Size
EV6100-4S1	SET consist of 1 (one) pouch, one o-ring, one belt and one drain tail closure.	Small ($\frac{7}{8}$ " Extra Tall	X-Small (300ml)
EV6110-4S1			X-Small (300ml)
EV6110-0S1		Medium ($\frac{7}{8}$ " - $1\frac{1}{8}$ " Extra Tall	Small (500ml)

STARTER Sets - Right Stoma Location - *Extra Large*

Product Number	Description	O-Ring Size	Pouch Size
EV6120-4S4	SET consist of 1 (one) pouch, one o-ring, one belt and one drain tail closure.	Extra Large ($1\frac{1}{4}$ "	X-Small (300ml)
EV6120-0S4			Small (500ml)

STANDARD Sets - Right Stoma Location

Product Number	Description	O-Ring Size	Pouch Size
EV-6200-004	SET consist of 2 (two) pouches, two o-rings, two belts and one drain tail closure.	Small ($\frac{7}{8}$ "	X-Small (300ml)
EV6200-000			Small (500ml)
EV6210-004		Medium ($\frac{7}{8}$ " - $1\frac{1}{8}$ "	X-Small (300ml)
EV6210-000			Small (500ml)
EV6220-004		Large ($1\frac{1}{8}$ " - $1\frac{3}{8}$ "	X-Small (300ml)
EV6220-000			Small (500ml)

Ileostomy Systems

STANDARD Sets - Right Stoma Location - Extra Tall

Product Number	Description	O-Ring Size	Pouch Size
EV6200-4S1	SETS consist of 2 (two) pouches, two o-rings, two belts and one drain tail closure.	Small ($\frac{7}{8}$ " Extra Tall	X-Small (300ml)
EV6210-0S1		Medium ($\frac{7}{8}$ " - 1- $\frac{1}{8}$ " Extra Tall	Small (500ml)
EV6220-4S1		Large (1 $\frac{1}{8}$ " - 1 $\frac{3}{8}$ " Extra Tall	X-Small (300ml)
EV6220-0S1			Small (500ml)

STANDARD Sets - Right Stoma Location - Extra Large

Product Number	Description	O-Ring Size	Pouch Size
EV6220-4S4	SETS consist of 2 (two) pouches, two o-rings, two belts and one drain tail closure.	Extra Large (1 $\frac{1}{4}$ "	X-Small (300ml)
EV6220-0S4			Small (500ml)

CONVENIENCE Sets - Right Stoma Location

Product Number	Description	O-Ring Size	Pouch Size
EV6400-004	SETS consist of 4 (four) pouches, two o-rings, two belts and one drain tail closure.	Small ($\frac{7}{8}$ "	X-Small (300ml)
EV6400-000			Small (500ml)
EV6410-004		Medium ($\frac{7}{8}$ " - 1- $\frac{1}{8}$ "	X-Small (300ml)
EV6410-000			Small (500ml)
EV6420-004		Large (1 $\frac{1}{8}$ " - 1 $\frac{3}{8}$ "	X-Small (300ml)
EV6420-000			Small (500ml)

CONVENIENCE Sets - Right Stoma Location - Extra Tall

Product Number	Description	O-Ring Size	Pouch Size
EV6400-4S1	SETS consist of 4 (four) pouches, two o-rings, two belts and one drain tail closure.	Small ($\frac{7}{8}$ " Extra Tall	X-Small (300ml)
EV6400-0S1			Small (500ml)
EV6410-4S1		Medium ($\frac{7}{8}$ " - 1- $\frac{1}{8}$ " Extra Tall	X-Small (300ml)
EV6410-0S1			Small (500ml)
EV6420-0S2		Large (1 $\frac{1}{8}$ " - 1 $\frac{3}{8}$ " Extra Tall	X-Small (300ml)
EV6420-0S1			Small (500ml)

Ileostomy Systems

CONVENIENCE Sets - Right Stoma Location - *Special Belt*

Product Number	Description	O-Ring Size	Pouch Size
EV6410-4S6	SETS consist of 4 (four) pouches, two o-rings, two special belts and one drain tail closure.	Medium ($\frac{7}{8}$ " - 1- $\frac{1}{8}$ ")	X-Small (300ml)

CONVENIENCE Sets - Right Stoma Location - *Extra Large*

Product Number	Description	O-Ring Size	Pouch Size
EV6420-0S4	SETS consist of 4 (four) pouches, two o-rings, two belts and one drain tail closure.	Extra Large (1 $\frac{1}{4}$ ")	Small (500ml)

STARTER Sets - Left Stoma Location

Product Number	Description	O-Ring Size	Pouch Size
EV6102-004	SETS consist of 1 (one) pouch, one o-ring, one belt and one drain tail closure.	Small ($\frac{7}{8}$ ")	X-Small (300ml)
EV6102-000			Small (500ml)
EV6112-004		Medium ($\frac{7}{8}$ " - 1- $\frac{1}{8}$ ")	X-Small (300ml)
EV6112-000			Small (500ml)
EV6122-004		Large (1 $\frac{1}{8}$ " - 1 $\frac{3}{8}$ ")	X-Small (300ml)
EV6122-000			Small (500ml)

STARTER Sets - Left Stoma Location - *Extra Tall*

Product Number	Description	O-Ring Size	Pouch Size
EV6102-4S1	SETS consist of 1 (one) pouch, one o-ring, one belt and one drain tail closure.	Small ($\frac{7}{8}$ ") Extra Tall	X-Small (300ml)
EV6112-0S1		Medium ($\frac{7}{8}$ " - 1- $\frac{1}{8}$ ") Extra Tall	Small (500ml)

STARTER Sets - Left Stoma Location - *Extra Large*

Product Number	Description	O-Ring Size	Pouch Size
EV6122-0S4	SETS consist of 1 (one) pouch, one o-ring, one belt and one drain tail closure.	Extra Large (1 $\frac{1}{4}$ ")	Small (500ml)

Ileostomy Systems

STANDARD Sets - Left Stoma Location

Product Number	Description	O-Ring Size	Pouch Size
EV6202-004	SETS consist of 2 (two) two pouches, two o-rings, two belts and one drain tail closure.	Small ($\frac{7}{8}$ "	X-Small (300ml)
EV6202-000			Small (500ml)
EV6212-004		Medium ($\frac{7}{8}$ " - $1\frac{1}{8}$ "	X-Small (300ml)
EV6212-000			Small (500ml)
EV6222-004		Large ($1\frac{1}{8}$ " - $1\frac{3}{8}$ "	X-Small (300ml)
EV6222-000			Small (500ml)

STANDARD Sets - Left Stoma Location - *Extra Tall*

Product Number	Description	O-Ring Size	Pouch Size
EV6212-0S1	SETS consist of 2 (two) pouches, two o-rings, two belts and one drain tail closure.	Medium ($\frac{7}{8}$ " - $1\frac{1}{8}$ ") Extra Tall	Small (500ml)
EV6222-0S1		Large ($1\frac{1}{8}$ " - $1\frac{3}{8}$ ") Extra Tall	Small (500ml)

STANDARD Sets - Left Stoma Location - *Extra Large*

Product Number	Description	O-Ring Size	Pouch Size
EV6222-0S4	SETS consist of 2 (two) pouches, two o-rings, two belts and one drain tail closure.	Extra Large ($1\frac{1}{4}$ "	Small (500ml)

CONVENIENCE Sets - Left Stoma Location

Product Number	Description	O-Ring Size	Pouch Size
EV6402-004	SETS consist of 4 (four) pouches, two o-rings, two belts and one drain tail closure.	Small ($\frac{7}{8}$ "	X-Small (300ml)
EV6402-000			Small (500ml)
EV6412-004		Medium ($\frac{7}{8}$ " - $1\frac{1}{8}$ "	X-Small (300ml)
EV6412-000			Small (500ml)
EV6422-004		Large ($1\frac{1}{8}$ " - $1\frac{3}{8}$ "	X-Small (300ml)
EV6422-000			Small (500ml)

Ileostomy Systems

CONVENIENCE Sets - Left Stoma Location - Extra Tall

Product Number	Description	O-Ring Size	Pouch Size
EV6412-OS1	SETS consist of 4 (four) pouches, two o-rings, two belts and one drain tail closure.	Medium ($\frac{7}{8}$ " - $1\frac{1}{8}$ ") Extra Tall	Small (500ml)
EV6422-OS1		Large ($1\frac{1}{8}$ " - $1\frac{3}{8}$ ") Extra Tall	Small (500ml)

CONVENIENCE Sets - Left Stoma Location - Extra Large



Product Number	Description	O-Ring Size	Pouch Size
EV6422-OS4	SETS consist of 4 (four) pouches, two o-rings, two belts and one drain tail closure.	Extra Large ($1\frac{1}{4}$ ")	Small (500ml)

Colostomy Systems

Open-Ended



Close-Ended



Colostomy Systems

STARTER Sets - Open-Ended

Product Number	Description	O-Ring Size	Pouch Flange Size
EV8250-0S3	SETS consist of 2 (two) pouches, one o-ring, one belt and one drain tail closure.	X-Small (1¼")	Small
EV8250-000		Small (1½")	Small
EV8260-000		Medium (2")	Medium
EV8270-000		Large (2½")	Large

STARTER Sets - Closed-Ended

Product Number	Description	O-Ring Size	Pouch Flange Size
EV8250-1S3	SETS consist of 2 (two) pouches, one o-ring, and one belt.	X-Small (1¼")	Small
EV8250-001		Small (1½")	Small
EV8260-001		Medium (2")	Medium
EV8270-001		Large (2½")	Large

CONVENIENCE Sets - Open-Ended

Product Number	Description	O-Ring Size	Pouch Flange Size
EV8050-0S3	SETS consist of 10 (ten) pouches, one o-ring, one belt and one drain tail closure.	X-Small (1¼")	Small
EV8050-000		Small (1½")	Small
EV8060-000		Medium (2")	Medium
EV8070-000		Large (2½")	Large

CONVENIENCE Sets - Closed-Ended

Product Number	Description	O-Ring Size	Pouch Flange Size
EV8050-1S3	SETS consist of 10 (ten) pouches, one o-ring, and one belt	X-Small (1¼")	Small
EV8050-001		Small (1½")	Small
EV8060-001		Medium (2")	Medium
EV8070-001		Large (2½")	Large



Silicone O-Rings




The pliable silicone o-ring, when tensioned properly by the belt, creates a secure seal between the pouch and abdomen.

*When added compression is needed to maintain a seal with a softer abdomen or skin roll situations, the Extra Tall ring may afford a more secure seal around the stoma.

Urostomy & Ileostomy

Product Number	Description	O-Ring Size	Fits Stoma Size
EV7000-00XX-RNG	Silicone O-Ring 1 (one) Per Package  Regular O-Ring (Side View)	Small	7/8 inch
EV7010-00XX-RNG		Medium	7/8 to 1-1/8 inch
EV7020-00XX-RNG		Large	1-1/8 - 1-3/8 inch
EV7020-01XX-RNG		X-Large	1-1/4 inch
EV7000-10XT-RNG	 *Extra-Tall O-Ring (Side View)	Small - Extra Tall*	7/8 inch
EV7010-10XT-RNG		Medium - Extra Tall*	7/8 to 1-1/8 inch
EV7020-10XT-RNG		Large - Extra Tall*	1-1/8 - 1-3/8 inch

Colostomy

Product Number	Description	O-Ring Size	Fits Stoma Size
EV7040-10XX-RNG	Silicone O-Ring 1 (one) Per Package  O-Ring (Side View)	X-Small	1-1/4 inch
EV7050-00XX-RNG		Small	1-1/2 inch
EV7060-00XX-RNG		Medium	2 inch
EV7070-00XX-RNG		Large	2-1/2 inch

Drain Tail Closures - Ileostomy & Colostomy



Product Number	Description
EV6030-00XX-CLS	Ileostomy & Colostomy Tail Closure

Belts



Belts are sized to match the pouch and stoma location.
So please order accordingly.

Belts are adjustable to fit waist sizes up to 55".
Inquire about larger sizes.

Urostomy

Product Number	Description	Fits Pouch Size	Stoma Location
EV7030-10XX-BLT	Belt holds urostomy product securely 1 (one) Per Package	X-Small (300ml)	Right & Left
EV7030-01XX-BLT		Small (500ml)	Right & Left
EV7030-04XX-BLT		Medium (800ml)	Right & Left
EV7030-07XX-BLT		Large (1100ml)	Right & Left


Ileostomy

Product Number	Description	Fits Pouch Size	Stoma Location
EV7030-18XX-BLT	Belt holds ileostomy product securely 1 (one) Per Package	X-Small (300ml)	Right & Left
EV7030-15XX-BLT		Small (500ml)	Right & Left

Colostomy

Product Number	Description	Belt Size
EV7030-25XX-BLT	Belt holds colostomy product securely 1 (one) Per Package	One size fits all, Right & Left

Pouch Comfort Panel (formally called Shields)



The soft cotton/polyester panel is used to provide a breathable layer between the pouch and body for added comfort.

Urostomy

Product Number	Description	Belt Size	Location
EV5033-004	Cotton/Polyester cloth Urostomy pouch comfort panel. 1 (one) Per Package	X-Small	Right & Left
EV5033-005		Small	Right & Left
EV5033-006		Medium	Right & Left
EV5033-007		Large	Right & Left

Ileostomy

Product Number	Description	Belt Size	Location
EV6033-04R	Cotton/Polyester cloth Ileostomy pouch comfort panel. 1 (one) Per Package	X-Small	Right
EV6033-04L		X-Small	Left
EV6033-05R		Small	Right
EV6033-05L		Small	Left

NU-HOPE LABORATORIES, INC
12640 Branford St., Pacoima, CA USA
(800) 899-5017
www.nu-hope.com

NHL1086